

Motivation

- Prostate cancer is the most common cancer affecting men.
- External beam radiation therapy is frequently used for treatment.
- Despite careful positioning at the start of treatment, the prostate may move during the dose delivery.

• Thus, there is a need for intrafraction monitoring of prostate position.

Non Rectilinear Frame

Images captured by the mechanical probe fall naturally in a cylindrical coordinate frame.



Probe: Side View









Volume: Side View

Top View





Conversion from Pixel to World Coordinates

Conversion from World to Pixel Coordinates

One limitation of the method is that for the cylindrical coordinate frame to remain valid, the probe must not move during acquisition.



Intrafraction prostate motion correction using a non-rectilinear image frame

Rupert Brooks, Elekta Ltd. Montreal, Quebec

TPUS prostate monitoring

We are interested in monitoring prostate position during treatment using 3D ultrasound taken trans-perineally (TPUS)





• TPUS imaging provides a non-invasive, non-ionizing means of prostate observation during radiation dose delivery.

Interpolation / Derivatives

Use of this coordinate frame in an intensity based registration framework requires the ability to interpolate image values at noninteger pixel positions, and to compute image derivatives relative to world space. $\frac{\partial \mathbf{I}}{\partial \mathbf{X}_{\mathbf{w}}} = \frac{\partial \mathbf{I}}{\partial \mathbf{X}_{\mathbf{p}}} \cdot \frac{\partial \mathbf{X}_{\mathbf{p}}}{\partial \mathbf{X}_{\mathbf{w}}}$





Interpolation is performed in the image space, which maps to a sector in physical space.



Image derivative $\left(\frac{\partial \mathbf{I}}{\partial \mathbf{X}_{w}}\right)$ with respect to world coordinates

$X_{p} = [i j$	k]: pixel coordinates	α:	ar
X _w :	world coordinates	r:	ra
X _o :	origin	s _i ,s _i :	pi
D:	matrix of direction cosines	u,v,w:	in

ngular step adial offset ixel scales in i and j directions ntermediate values

Selected Results on Typical Series

rupert.brooks@elekta.com



Equipment

• 3D ultrasound sweeps were obtained using a m4DC7-3/40 Convex 4D probe (Ultrasonix, Vancouver, Canada) mechanically swept probe.



• Tracked images were captured and stored for offline processing using a modified version of the Elekta Clarity US-Guide system

Evaluation and Discussion

Image sequences were captured from healthy volunteers. The prostate position was tracked using sequential registration. The direct use of a non-rectilinear image was compared to results obtained with a typical reconstruct and register approach.

		Reconstructed		Cylindrical		Centroid	
Case	# Vols	Time	σ	, Time	σ	mTRE	
P1S1	115	0.277	0.143	0.827	0.401	0.160	
P1S2	115	0.197	0.102	0.533	0.273	0.636	
P2S2	115	0.222	0.075	0.638	0.284	0.193	
P2S3	115	0.208	0.069	0.562	0.206	0.265	
P2S4	131	0.193	0.053	0.551	0.146	0.113	
P3S1	115	0.306	0.203	0.864	0.424	0.619	
Overall	706	0.233	0.118	0.660	0.303	0.320	

Mean registration time and positional difference between methods. **Reconstruction requires an additional 2 seconds per image**

These results show that registration using cylindrical coordinates compares favorably with the use of a rectilinear image. While registration with the cylindrical image is slightly slower, there is a considerable overall time saving by avoiding reconstruction.













